



QUALITY FOR LIFE

Start Junior

Lightweight Wheelchair for Children and Adolescents

Maximum User Weight 90kg
Article Number 480F53=60000

Company Name _____

Contact Person _____

Customer Name _____

Delivery To _____

Purchase Order No _____



Standard Equipment

- Aluminium Frame
- Integrated Tip Assist
- Individual black swing in/out, angle adjustable plastic footplates
- Flip up, removable armrest- 300mm pads, adjustable position with fixed height adjustment
- Rear wheel with quick release axle
- Standard wheel lock
- Standard Nylon back upholstery
- Lower leg length from 16cm
- Seat depth adjustable from 31 to 37cm
- Rear seat to floor height from 37cm to 52cm (see Instructions for Use for matrix)
- Back Height adjustable from 30cm to 37.5cm

Selection	[]	[]	[]
Seat Width	28	30.5	35.5
Axle plate	Standard	Standard	Standard
Frame Coating	Red	Red	Red
Seat Depth	34 (adjustable from 31 to 37cm)	34 (adjustable from 31 to 37cm)	34 (adjustable from 31 to 37cm)
Footrest	Short (16-30cm)	Short (16-30cm)	Standard (31cm+)

Otto Bock Australia Pty. Ltd.

ABN 35 001 242 360

PO Box 7224 Baulkham Hills Business Centre N.S.W. 2153 AUSTRALIA

Suite 1.01 Century Corporate Centre, 62 Norwest Boulevard, Norwest Business Park, Baulkham Hills N.S.W. 2153 AUSTRALIA

Phone: +61 2 8818 2800 • Fax: +61 2 8814 4500 • Email: healthcare@ottobock.com.au

www.ottobock.com.au



QUALITY FOR LIFE

Back Height	32.5 (adjustable from 30 to 37.5cm)	32.5 (adjustable from 30 to 37.5cm)	32.5 (adjustable from 30 to 37.5cm)
Back Tube	Lumbar Curve	Lumbar curve	Lumbar curve
Push Handles	Standard	Standard	Standard
Back Fabric	Standard	Standard	Standard
Caster Size	5"	5"	6"
Rear Wheel	20" PU	22" PU	22" PU
Side Panel	Short Mounted	Short Mounted	Short Mounted
Armrest	standard	standard	standard
Caster Fork	standard	standard	standard

Additional Accessory Options

<input type="checkbox"/>	Height adjustable, angled push handles		481D53=PK050
<input type="checkbox"/>	Crutch Holder	<input type="checkbox"/> right <input type="checkbox"/> left	481A53=SK018
<input type="checkbox"/>	Anti Tippers (each)		481A53=SK070
<input type="checkbox"/>	Elevating Legrests	<input type="checkbox"/> right <input type="checkbox"/> left	481B53=RK126 481B53=LK126
<input type="checkbox"/>	Heel Bands		481B53=SK010
<input type="checkbox"/>	Plug on Wheel lock lever extension (each)		481H53=SK012
<input type="checkbox"/>	Cubic Foam cushion	(depth=___ cm)	476C00=SK100
<input type="checkbox"/>	Combi Foam cushion	(depth=___ cm)	476C00=SK200
<input type="checkbox"/>	EvoLight cushion	(depth=___ cm)	476C00=SK350
<input type="checkbox"/>	Pelvic belt		446B1=2-7
<input type="checkbox"/>	22" pneumatic tyre and tube (unfitted)		TBA
<input type="checkbox"/>	24" pneumatic tyre and tube (unfitted)		481G09=SK010-WH
<input type="checkbox"/>	Therapy tray (slide on)		481E53=SK0**
<input type="checkbox"/>	Spoke Protector (specify design_____)		481G00=SK7**
<input type="checkbox"/>	Rear wheel with drum brake(not available for 20" wheels)		481G53=PK0**
<input type="checkbox"/>	One arm drive (not available for 20" wheels)		481G53=SK0**

Otto Bock Australia Pty. Ltd.

ABN 35 001 242 360

PO Box 7224 Baulkham Hills Business Centre N.S.W. 2153 AUSTRALIA

Suite 1.01 Century Corporate Centre, 62 Norwest Boulevard, Norwest Business Park, Baulkham Hills N.S.W. 2153 AUSTRALIA

Phone: +61 2 8818 2800 • Fax: +61 2 8814 4500 • Email: healthcare@ottobock.com.au

www.ottobock.com.au